

Purified anti-SSRP1

Mouse IgG1, kappa
3E4



BioLegend
From Revvity

100 µg

Artikel-Nr.: B609802 | **BioLegend** | Hersteller-Nr.: 609802

175,00 €*

*zzgl. MwSt. zzgl. Versandkosten

Beschreibung

App_Abbrev: FA, FC, ICC, ICFC, IHC, IP, WB
Clone: 3E4
Format: PURE
Isotype: Mouse IgG1, kappa
Species of Origin: Mouse
Verpackung: 100 µg

Purified anti-SSRP1

Clone: 3E4

Isotype: Mouse IgG1, kappa

Other_Names: Structure specific recognition protein 1 (SSRP1), Facilitates Chromatin Transcription 80 kD subunit, FACTp80, Recombination signal sequence recognition protein

Size: 100 µg

Reactivity: Human

Format: PURE

App_Abbrev: FC, ICFC, ICC, IHC, IP, WB, FA

Species of Origin: Mouse

Weitere Informationen zum Produkt. Link zum Hersteller:

<https://www.biolegend.com/Default.aspx?ID=6664&productid=1991>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>