

Purified anti-Spt16 (FACT140 complex)

Mouse IgG2a, kappa
8D2



BioLegend
From Revvity

25 µg

Artikel-Nr.: B607001 | [BioLegend](#) | Hersteller-Nr.: 607001

81,00 €*

*zzgl. MwSt. [zzgl. Versandkosten](#)

Beschreibung

App_Abbrev: FA, FC, ICC, ICFC, IHC, IP, WB
Clone: 8D2
Format: PURE
Isotype: Mouse IgG2a, kappa
Species of Origin: Mouse
Verpackung: 25 µg

Purified anti-Spt16 (FACT140 complex)

Clone: 8D2

Isotype: Mouse IgG2a, kappa

Other_Names: Facilitates Chromatin Transcription (FACT140), Spt 16 protein, Cdc68

Size: 25 µg

Reactivity: Human, Mouse

Format: PURE

App_Abbrev: FC, ICFC, ICC, IHC, IP, WB, FA

Species of Origin: Mouse

Weitere Informationen zum Produkt. Link zum Hersteller:

<https://www.biolegend.com/Default.aspx?ID=6664&productid=2003>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>