

Purified anti-SSRP1

Mouse IgG2b, kappa
10D1



BioLegend
From Revvity

25 µg
Artikel-Nr.: B609701 | [BioLegend](#) | Hersteller-Nr.: 609701

81,00 €*

*zzgl. MwSt. [zzgl. Versandkosten](#)

Beschreibung

App_Abbrev: FA, FC, ICC, ICFC, IHC, IP, WB
Clone: 10D1
Format: PURE
Isotype: Mouse IgG2b, kappa
Species of Origin: Mouse
Verpackung: 25 µg

Purified anti-SSRP1

Clone: 10D1

Isotype: Mouse IgG2b, kappa

Other_Names: Structure specific recognition protein 1 (SSRP1), Facilitates Chromatin Transcription 80 kD subunit, FACTp80, Recombination signal sequence recognition protein

Size: 25 µg

Reactivity: Human, Mouse

Format: PURE

App_Abbrev: FC, ICFC, ICC, IHC, IP, WB, FA

Species of Origin: Mouse

Weitere Informationen zum Produkt. Link zum Hersteller:

<https://www.biolegend.com/Default.aspx?ID=6664&productid=1978>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>