

Purified anti-C/EBP beta

Mouse IgG2a, kappa
1H7



BioLegend
From Revvity

100 µg
Artikel-Nr.: B606202 | [BioLegend](#) | Hersteller-Nr.: 606202

175,00 €*

*zzgl. MwSt. [zzgl. Versandkosten](#)

Beschreibung

App_Abbrev: FA, FC, ICC, ICFC, IHC, IP, WB
Clone: 1H7
Format: PURE
Isotype: Mouse IgG2a, kappa
Species of Origin: Mouse
Verpackung: 100 µg

Purified anti-C/EBP beta

Clone: 1H7

Isotype: Mouse IgG2a, kappa

Other_Names: CCAAT/enhancer-binding protein beta, Alpha-1-acid glycoprotein/enhancer binding protein, Interleukin 6-dependent DNA-binding protein (IL-6DBP), Nuclear factor of interleukin 6 (NF-IL6), Liver-enriched transcriptional activator protein (LAP)

Size: 100 µg

Reactivity: Mouse, Rat

Format: PURE

App_Abbrev: FC, ICFC, ICC, IHC, IP, WB, FA

Species of Origin: Mouse

Weitere Informationen zum Produkt. Link zum Hersteller:

<https://www.biolegend.com/Default.aspx?ID=6664&productid=1984>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>