

Biotin anti-human IL-1 RA

Rabbit Polyclonal IgG
Poly5095



BioLegend
From Revvity

500 µg
Artikel-Nr.: B509502 | **BioLegend** | Hersteller-Nr.: 509502

300,00 €*

*zzgl. MwSt. zzgl. Versandkosten

Beschreibung

App_Abbrev: FC, ICFC
Clone: Poly5095
Format: BIOTIN
Isotype: Rabbit Polyclonal IgG
Species of Origin: Rabbit
Verpackung: 500 µg

Biotin anti-human IL-1 RA

Clone: Poly5095

Isotype: Rabbit Polyclonal IgG

Other_Names: Interleukin-1 receptor antagonist, IL-1 inhibitor, IL-1 receptor antagonist protein (IRAP), IL-1 inhibitory activity, ILS, MG50, IL-1RN

Size: 500 µg

Reactivity: Human

Format: BIOTIN

App_Abbrev: FC, ICFC

Species of Origin: Rabbit

Weitere Informationen zum Produkt. Link zum Hersteller:

<https://www.biolegend.com/Default.aspx?ID=6664&productid=1789>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>