

Purified anti-human IL-1 RA

Mouse IgG1
JK1RA-1



BioLegend
From Revvity

50 µg
Artikel-Nr.: B509901 | **BioLegend** | Hersteller-Nr.: 509901

85,00 €*

*zzgl. MwSt. zzgl. Versandkosten

Beschreibung

App_Abbrev: FA, FC, ICC, ICFC, IHC, IP, WB
Clone: JK1RA-1
Format: PURE
Isotype: Mouse IgG1
Species of Origin: Mouse
Verpackung: 50 µg

Purified anti-human IL-1 RA

Clone: JK1RA-1

Isotype: Mouse IgG1

Other_Names: Interleukin-1 receptor antagonist, IL-1 inhibitor, IL-1 receptor antagonist protein (IRAP), IL-1 inhibitory activity, ILS, MG50, IL-1RN

Size: 50 µg

Reactivity: Human

Format: PURE

App_Abbrev: FC, ICFC, ICC, IHC, IP, WB, FA

Species of Origin: Mouse

Weitere Informationen zum Produkt. Link zum Hersteller:

<https://www.biolegend.com/Default.aspx?ID=6664&productid=1788>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>