

Purified anti-mouse CD254 (TRANCE, RANKL)

Rat IgG2a, kappa
IK22/5



BioLegend
From Revvity

500 µg
Artikel-Nr.: B510002 | **BioLegend** | Hersteller-Nr.: 510002

109,00 €*

*zzgl. MwSt. zzgl. Versandkosten

Beschreibung

App_Abbrev: FA, FC, ICC, ICFC, IHC, IP, WB
Clone: IK22/5
Format: PURE
Isotype: Rat IgG2a, kappa
Species of Origin: Rat
Verpackung: 500 µg

Purified anti-mouse CD254 (TRANCE, RANKL)

Clone: IK22/5

Isotype: Rat IgG2a, kappa

Other_Names: Tumor necrosis factor-related activation-induced cytokine, Receptor activator of NF-kappaB ligand (RANK L), Osteoprotegerin ligand (OPGL), Osteoclast differentiation factor (ODF), SOFA, TNFSF-11, CD254, RANKL

Size: 500 µg

Reactivity: Mouse

Format: PURE

App_Abbrev: FC, ICFC, ICC, IHC, IP, WB, FA

Species of Origin: Rat

Weitere Informationen zum Produkt. Link zum Hersteller: <https://www.biolegend.com/Default.aspx?ID=6664&productid=310>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>