

PE anti-mouse CD254 (TRANCE, RANKL)

Rat IgG2a, kappa
IK22/5



BioLegend
From Revvity

50 µg
Artikel-Nr.: B510005 | **BioLegend** | Hersteller-Nr.: 510005

156,00 €*

*zzgl. MwSt. zzgl. Versandkosten

Beschreibung

App_Abbrev: FC, ICFC
Clone: IK22/5
Format: PE
Isotype: Rat IgG2a, kappa
Species of Origin: Rat
Verpackung: 50 µg

PE anti-mouse CD254 (TRANCE, RANKL)

Clone: IK22/5

Isotype: Rat IgG2a, kappa

Other_Names: Tumor necrosis factor-related activation-induced cytokine, Receptor activator of NF-kappaB ligand (RANK L), Osteoprotegerin ligand (OPGL), Osteoclast differentiation factor (ODF), SOFA, TNFSF-11, CD254, RANKL

Size: 50 µg

Reactivity: Mouse

Format: PE

App_Abbrev: FC, ICFC

Species of Origin: Rat

Weitere Informationen zum Produkt. Link zum Hersteller: <https://www.biolegend.com/Default.aspx?ID=6664&productid=309>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>