

Purified anti-human IL-23 (p19)

Mouse IgG1, kappa
HLT2736



BioLegend
From Revvity

50 µg
Artikel-Nr.: B511201 | [BioLegend](#) | Hersteller-Nr.: 511201

95,00 €*

*zzgl. MwSt. [zzgl. Versandkosten](#)

Beschreibung

App_Abbrev: FA, FC, ICC, ICFC, IHC, IP, WB
Clone: HLT2736
Format: PURE
Isotype: Mouse IgG1, kappa
Species of Origin: Mouse
Verpackung: 50 µg

Purified anti-human IL-23 (p19)

Clone: HLT2736
Isotype: Mouse IgG1, kappa
Other_Names: p19p40, JKA3 induced upon T cell activation
Size: 50 µg

Reactivity: Human
Format: PURE
App_Abbrev: FC, ICFC, ICC, IHC, IP, WB, FA
Species of Origin: Mouse

Weitere Informationen zum Produkt. Link zum Hersteller:
<https://www.biolegend.com/Default.aspx?ID=6664&productid=2842>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>