

Purified anti-human IL-9

Mouse IgG2b, kappa
MH9A4



BioLegend
From Revvity

500 µg

Artikel-Nr.: B507602 | **BioLegend** | Hersteller-Nr.: 507602

226,00 €*

*zzgl. MwSt. zzgl. Versandkosten

Beschreibung

App_Abbrev: FA, FC, ICC, ICFC, IHC, IP, WB
Clone: MH9A4
Format: PURE
Isotype: Mouse IgG2b, kappa
Species of Origin: Mouse
Verpackung: 500 µg

Purified anti-human IL-9

Clone: MH9A4

Isotype: Mouse IgG2b, kappa

Other_Names: Interleukin-9, P40, Mast cell growth factor (MCGF), Megakaryoblast growth factor, T cell growth factor-3 (TCGF-3)

Size: 500 µg

Reactivity: Human

Format: PURE

App_Abbrev: FC, ICFC, ICC, IHC, IP, WB, FA

Species of Origin: Mouse

Weitere Informationen zum Produkt. Link zum Hersteller:

<https://www.biolegend.com/Default.aspx?ID=6664&productid=1621>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>