

PE anti-human IL-9

Mouse IgG2b, kappa
MH9A4



BioLegend
From Revvity

50 µg
Artikel-Nr.: B507603 | **BioLegend** | Hersteller-Nr.: 507603

147,00 €*

*zzgl. MwSt. zzgl. Versandkosten

Beschreibung

App_Abbrev: FC, ICFC
Clone: MH9A4
Format: PE
Isotype: Mouse IgG2b, kappa
Species of Origin: Mouse
Verpackung: 50 µg

PE anti-human IL-9

Clone: MH9A4

Isotype: Mouse IgG2b, kappa

Other_Names: Interleukin-9, P40, Mast cell growth factor (MCGF), Megakaryoblast growth factor, T cell growth factor-3 (TCGF-3)

Size: 50 µg

Reactivity: Human

Format: PE

App_Abbrev: FC, ICFC

Species of Origin: Mouse

Weitere Informationen zum Produkt. Link zum Hersteller:

<https://www.biolegend.com/Default.aspx?ID=6664&productid=1974>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>