

Purified anti-mouse IFN-beta

Armenian Hamster IgG
MIB-5E9.1



BioLegend
From Revvity

500 µg
Artikel-Nr.: B508102 | **BioLegend** | Hersteller-Nr.: 508102

226,00 €*

*zzgl. MwSt. zzgl. Versandkosten

Beschreibung

App_Abbrev: FA, FC, ICC, ICFC, IHC, IP, WB
Clone: MIB-5E9.1
Format: PURE
Isotype: Armenian Hamster IgG
Species of Origin: Armenian Hamster
Verpackung: 500 µg

Purified anti-mouse IFN-beta

Clone: MIB-5E9.1

Isotype: Armenian Hamster IgG

Other_Names: Interferon-beta, Fibroblast interferon (Fi-IFN, F-IFN), Type-1 interferon, acid-stable interferon, R1-G1 factor, IFN-b, IFN-beta

Size: 500 µg

Reactivity: Mouse

Format: PURE

App_Abbrev: FC, ICFC, ICC, IHC, IP, WB, FA

Species of Origin: Armenian Hamster

Weitere Informationen zum Produkt. Link zum Hersteller:

<https://www.biolegend.com/Default.aspx?ID=6664&productid=1623>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>