

Purified anti-human IL-1beta

Mouse IgG2b, kappa
JK1B-2



BioLegend
From Revvity

500 µg
Artikel-Nr.: B508304 | **BioLegend** | Hersteller-Nr.: 508304

278,00 €*

*zzgl. MwSt. zzgl. Versandkosten

Beschreibung

App_Abbrev: FA, FC, ICC, ICFC, IHC, IP, WB
Clone: JK1B-2
Format: PURE
Isotype: Mouse IgG2b, kappa
Species of Origin: Mouse
Verpackung: 500 µg

Purified anti-human IL-1beta

Clone: JK1B-2

Isotype: Mouse IgG2b, kappa

Other_Names: Interleukin-1beta, Catabolin, Hematopoietin-1 (H1), IFNbeta-inducing factor, Interleukin-beta, Osteoclast activating factor (OAF)

Size: 500 µg

Reactivity: Human

Format: PURE

App_Abbrev: FC, ICFC, ICC, IHC, IP, WB, FA

Species of Origin: Mouse

Weitere Informationen zum Produkt. Link zum Hersteller:

<https://www.biolegend.com/Default.aspx?ID=6664&productid=2493>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>