

Biotin anti-human IL-12/IL-23 p40

Mouse IgG1, kappa
C8.6



BioLegend
From Revvity

500 µg
Artikel-Nr.: B508802 | **BioLegend** | Hersteller-Nr.: 508802

220,00 €*

*zzgl. MwSt. zzgl. Versandkosten

Beschreibung

App_Abbrev: FC, ICFC
Clone: C8.6
Format: BIOTIN
Isotype: Mouse IgG1, kappa
Species of Origin: Mouse
Verpackung: 500 µg

Biotin anti-human IL-12/IL-23 p40

Clone: C8.6

Isotype: Mouse IgG1, kappa

Other_Names: Interleukin-12 p40, Interleukin-23 p40, Cytotoxic lymphocyte maturation factor (CLMF), Natural killer cell stimulatory factor (NKSF), CTL maturation factor (TcMF), T-cell stimulating factor (TSF)

Size: 500 µg

Reactivity: Human

Format: BIOTIN

App_Abbrev: FC, ICFC

Species of Origin: Mouse

Weitere Informationen zum Produkt. Link zum Hersteller:

<https://www.biolegend.com/Default.aspx?ID=6664&productid=1008>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>