

Alexa Fluor® 700 anti-human IFN-gamma

Mouse IgG1, kappa
B27



BioLegend
From Revvity

25 µg
Artikel-Nr.: B506515 | **BioLegend** | Hersteller-Nr.: 506515

101,00 €*

*zzgl. MwSt. zzgl. Versandkosten

Beschreibung

App_Abbrev: FC, ICFC
Clone: B27
Format: A700
Isotype: Mouse IgG1, kappa
Species of Origin: Mouse
Verpackung: 25 µg

Alexa Fluor® 700 anti-human IFN-gamma

Clone: B27

Isotype: Mouse IgG1, kappa

Other_Names: Interferon-gamma, Immune interferon, Type II interferon, T cell interferon, Macrophage-activating factor (MAF), IFN-g, IFN-gamma

Size: 25 µg

Reactivity: Human

Format: A700

App_Abbrev: FC, ICFC

Species of Origin: Mouse

Weitere Informationen zum Produkt. Link zum Hersteller:

<https://www.biolegend.com/Default.aspx?ID=6664&productid=3445>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>