

Purified anti-human IL-7

Rat IgG2a, kappa
BVD10-11C10



BioLegend
From Revvity

500 µg

Artikel-Nr.: B506604 | [BioLegend](#) | Hersteller-Nr.: 506604

235,00 €*

*zzgl. MwSt. [zzgl. Versandkosten](#)

Beschreibung

App_Abbrev: FA, FC, ICC, ICFC, IHC, IP, WB
Clone: BVD10-11C10
Format: PURE
Isotype: Rat IgG2a, kappa
Species of Origin: Rat
Verpackung: 500 µg

Purified anti-human IL-7

Clone: BVD10-11C10

Isotype: Rat IgG2a, kappa

Other_Names: Interleukin-7, Lymphopoietin-1 (LP-1), Pre-B cell growth factor, Thymocyte growth factor

Size: 500 µg

Reactivity: Human

Format: PURE

App_Abbrev: FC, ICFC, ICC, IHC, IP, WB, FA

Species of Origin: Rat

Weitere Informationen zum Produkt. Link zum Hersteller:

<https://www.biolegend.com/Default.aspx?ID=6664&productid=3680>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>