

## Purified anti-human IL-10

Rat IgG2a, kappa  
JES3-19F1



**BioLegend**  
From Revvity

50 µg  
Artikel-Nr.: B506801 | [BioLegend](#) | Hersteller-Nr.: 506801

**67,00 €\***

\*zzgl. MwSt. [zzgl. Versandkosten](#)

### Beschreibung

**App\_Abbrev:** FA, FC, ICC, ICFC, IHC, IP, WB  
**Clone:** JES3-19F1  
**Format:** PURE  
**Isotype:** Rat IgG2a, kappa  
**Species of Origin:** Rat  
**Verpackung:** 50 µg

Purified anti-human IL-10

Clone: JES3-19F1

Isotype: Rat IgG2a, kappa

Other\_Names: Interleukin 10, B cell derived T cell growth factor (B-TCGF), Cytokine synthesis inhibitory factor (CSIF), T-cell growth inhibitory factor (TGIF)

Size: 50 µg

Reactivity: Human

Format: PURE

App\_Abbrev: FC, ICFC, ICC, IHC, IP, WB, FA

Species of Origin: Rat

Weitere Informationen zum Produkt. Link zum Hersteller:

<https://www.biolegend.com/Default.aspx?ID=6664&productid=1575>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>