

APC anti-human IL-10

Rat IgG2a, kappa
JES3-19F1



BioLegend
From Revvity

100 tests

Artikel-Nr.: B506807 | [BioLegend](#) | Hersteller-Nr.: 506807

237,00 €*

*zzgl. MwSt. [zzgl. Versandkosten](#)

Beschreibung

App_Abbrev: FC, ICFC
Clone: JES3-19F1
Format: APC
Isotype: Rat IgG2a, kappa
Species of Origin: Rat
Verpackung: 100 Tests

APC anti-human IL-10

Clone: JES3-19F1

Isotype: Rat IgG2a, kappa

Other_Names: Interleukin 10, B cell derived T cell growth factor (B-TCGF), Cytokine synthesis inhibitory factor (CSIF), T-cell growth inhibitory factor (TGIF)

Size: 100 tests

Reactivity: Human

Format: APC

App_Abbrev: FC, ICFC

Species of Origin: Rat

Weitere Informationen zum Produkt. Link zum Hersteller:

<https://www.biolegend.com/Default.aspx?ID=6664&productid=1572>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>