

## PE anti-mouse IL-17A

Rat IgG1, kappa  
TC11-18H10.1



**BioLegend**  
From Revvity

25 µg  
Artikel-Nr.: B506903 | **BioLegend** | Hersteller-Nr.: 506903

**88,00 €\***

\*zzgl. MwSt. zzgl. Versandkosten

### Beschreibung

**App\_Abbrev:** FC, ICFC  
**Clone:** TC11-18H10.1  
**Format:** PE  
**Isotype:** Rat IgG1, kappa  
**Species of Origin:** Rat  
**Verpackung:** 25 µg

PE anti-mouse IL-17A

Clone: TC11-18H10.1  
Isotype: Rat IgG1, kappa  
Other\_Names: Interleukin-17, Cytotoxic T lymphocyte-associated antigen 8 (CTLA-8)  
Size: 25 µg

Reactivity: Mouse  
Format: PE  
App\_Abbrev: FC, ICFC  
Species of Origin: Rat

Weitere Informationen zum Produkt. Link zum Hersteller:  
<https://www.biolegend.com/Default.aspx?ID=6664&productid=1633>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>