

PE anti-mouse IL-17A

Rat IgG1, kappa
TC11-18H10.1



BioLegend
From Revvity

100 µg

Artikel-Nr.: B506904 | **BioLegend** | Hersteller-Nr.: 506904

224,00 €*

*zzgl. MwSt. zzgl. Versandkosten

Beschreibung

App_Abbrev: FC, ICFC
Clone: TC11-18H10.1
Format: PE
Isotype: Rat IgG1, kappa
Species of Origin: Rat
Verpackung: 100 µg

PE anti-mouse IL-17A

Clone: TC11-18H10.1

Isotype: Rat IgG1, kappa

Other_Names: Interleukin-17, Cytotoxic T lymphocyte-associated antigen 8 (CTLA-8)

Size: 100 µg

Reactivity: Mouse

Format: PE

App_Abbrev: FC, ICFC

Species of Origin: Rat

Weitere Informationen zum Produkt. Link zum Hersteller:

<https://www.biolegend.com/Default.aspx?ID=6664&productid=1633>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>