

APC anti-mouse IL-17A

Rat IgG1, kappa
TC11-18H10.1



BioLegend
From Revvity

25 µg
Artikel-Nr.: B506915 | **BioLegend** | Hersteller-Nr.: 506915

105,00 €*

*zzgl. MwSt. zzgl. Versandkosten

Beschreibung

App_Abbrev: FC, ICFC
Clone: TC11-18H10.1
Format: APC
Isotype: Rat IgG1, kappa
Species of Origin: Rat
Verpackung: 25 µg

APC anti-mouse IL-17A

Clone: TC11-18H10.1
Isotype: Rat IgG1, kappa
Other_Names: Interleukin-17, Cytotoxic T lymphocyte-associated antigen 8 (CTLA-8)
Size: 25 µg

Reactivity: Mouse
Format: APC
App_Abbrev: FC, ICFC
Species of Origin: Rat

Weitere Informationen zum Produkt. Link zum Hersteller:
<https://www.biolegend.com/Default.aspx?ID=6664&productid=3540>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>