

Purified anti-mouse GM-CSF

Rat IgG1, kappa
MP1-31G6



BioLegend
From Revvity

500 µg
Artikel-Nr.: B505504 | **BioLegend** | Hersteller-Nr.: 505504

188,00 €*

*zzgl. MwSt. zzgl. Versandkosten

Beschreibung

App_Abbrev: FA, FC, ICC, ICFC, IHC, IP, WB
Clone: MP1-31G6
Format: PURE
Isotype: Rat IgG1, kappa
Species of Origin: Rat
Verpackung: 500 µg

Purified anti-mouse GM-CSF

Clone: MP1-31G6

Isotype: Rat IgG1, kappa

Other_Names: Granulocyte/macrophage-colony stimulating factor, CSF-alpha, Pluripoietin-alpha, Eosinophil colony stimulating factor (Eo-CSF), Burst promoting activity (BPA)

Size: 500 µg

Reactivity: Mouse

Format: PURE

App_Abbrev: FC, ICFC, ICC, IHC, IP, WB, FA

Species of Origin: Rat

Weitere Informationen zum Produkt. Link zum Hersteller:

<https://www.biolegend.com/Default.aspx?ID=6664&productid=2220>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>