

Purified anti-mouse IFN-gamma

Rat IgG1, kappa
R4-6A2



BioLegend
From Revvity

50 µg

Artikel-Nr.: B505701 | **BioLegend** | Hersteller-Nr.: 505701

67,00 €*

*zzgl. MwSt. zzgl. Versandkosten

Beschreibung

App_Abbrev: FA, FC, ICC, ICFC, IHC, IP, WB
Clone: R4-6A2
Format: PURE
Isotype: Rat IgG1, kappa
Species of Origin: Rat
Verpackung: 50 µg

Purified anti-mouse IFN-gamma

Clone: R4-6A2

Isotype: Rat IgG1, kappa

Other_Names: Interferon-gamma, Immune interferon, Type II interferon, T cell interferon, Macrophage-activating factor (MAF), IFN-g, IFN-gamma

Size: 50 µg

Reactivity: Mouse

Format: PURE

App_Abbrev: FC, ICFC, ICC, IHC, IP, WB, FA

Species of Origin: Rat

Weitere Informationen zum Produkt. Link zum Hersteller: <https://www.biolegend.com/Default.aspx?ID=6664&productid=987>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>