

Alexa Fluor® 647 anti-mouse IFN-gamma

Rat IgG1, kappa
XMG1.2



BioLegend
From Revvity

25 µg

Artikel-Nr.: B505816 | [BioLegend](#) | Hersteller-Nr.: 505816

90,00 €*

*zzgl. MwSt. [zzgl. Versandkosten](#)

Beschreibung

App_Abbrev: FC, ICFC
Clone: XMG1.2
Format: A647
Isotype: Rat IgG1, kappa
Species of Origin: Rat
Verpackung: 25 µg

Alexa Fluor® 647 anti-mouse IFN-gamma

Clone: XMG1.2

Isotype: Rat IgG1, kappa

Other_Names: Interferon-gamma, Immune interferon, Type II interferon, T cell interferon, Macrophage-activating factor (MAF)

Size: 25 µg

Reactivity: Mouse

Format: A647

App_Abbrev: FC, ICFC

Species of Origin: Rat

Weitere Informationen zum Produkt. Link zum Hersteller:

<https://www.biolegend.com/Default.aspx?ID=6664&productid=2722>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>