

Purified anti-mouse/rat/human MCP-1

Hamster IgG
2H5



BioLegend
From Revvity

50 µg
Artikel-Nr.: B505901 | **BioLegend** | Hersteller-Nr.: 505901

76,00 €*

*zzgl. MwSt. zzgl. Versandkosten

Beschreibung

App_Abbrev: FA, FC, ICC, ICFC, IHC, IP, WB
Clone: 2H5
Format: PURE
Isotype: Hamster IgG
Species of Origin: Hamster
Verpackung: 50 µg

Purified anti-mouse/rat/human MCP-1

Clone: 2H5

Isotype: Hamster IgG

Other_Names: Monocyte chemoattractant protein-1, Monocyte chemoattractant and activating factor (MCAF), JE, Small inducible cytokine A2 (SCYA2), HC-11, P6, Smooth muscle cell chemotactic factor (SMC-CF), CCL2

Size: 50 µg

Reactivity: Mouse, Rat, Human

Format: PURE

App_Abbrev: FC, ICFC, ICC, IHC, IP, WB, FA

Species of Origin: Hamster

Weitere Informationen zum Produkt. Link zum Hersteller:

<https://www.biolegend.com/Default.aspx?ID=6664&productid=1003>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>