

Purified anti-mouse/rat TNF-alpha

Armenian Hamster IgG
TN3-19.12



BioLegend
From Revvity

50 µg
Artikel-Nr.: B506101 | **BioLegend** | Hersteller-Nr.: 506101

81,00 €*

*zzgl. MwSt. zzgl. Versandkosten

Beschreibung

App_Abbrev: FA, FC, ICC, ICFC, IHC, IP, WB
Clone: TN3-19.12
Format: PURE
Isotype: Armenian Hamster IgG
Species of Origin: Armenian Hamster
Verpackung: 50 µg

Purified anti-mouse/rat TNF-alpha

Clone: TN3-19.12

Isotype: Armenian Hamster IgG

Other_Names: Tumor necrosis factor-alpha, Cachectin, Necrosin, Macrophage cytotoxic factor (MCF), Differentiation inducing factor (DIF), TNFSF-2

Size: 50 µg

Reactivity: Mouse, Rat

Format: PURE

App_Abbrev: FC, ICFC, ICC, IHC, IP, WB, FA

Species of Origin: Armenian Hamster

Weitere Informationen zum Produkt. Link zum Hersteller:

<https://www.biolegend.com/Default.aspx?ID=6664&productid=1367>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>