

Alexa Fluor® 488 anti-mouse TNF-alpha

Rat IgG1, kappa
MP6-XT22



BioLegend
From Revvity

100 µg
Artikel-Nr.: B506313 | **BioLegend** | Hersteller-Nr.: 506313

215,00 €*

*zzgl. MwSt. zzgl. Versandkosten

Beschreibung

App_Abbrev: FC, ICFC
Clone: MP6-XT22
Format: A488
Isotype: Rat IgG1, kappa
Species of Origin: Rat
Verpackung: 100 µg

Alexa Fluor® 488 anti-mouse TNF-alpha

Clone: MP6-XT22

Isotype: Rat IgG1, kappa

Other_Names: Tumor necrosis factor-alpha, Cachectin, Necrosin, Macrophage cytotoxic factor (MCF), Differentiation inducing factor (DIF), TNFSF-2, TNF-a, TNF-alpha

Size: 100 µg

Reactivity: Mouse

Format: A488

App_Abbrev: FC, ICFC

Species of Origin: Rat

Weitere Informationen zum Produkt. Link zum Hersteller:

<https://www.biolegend.com/Default.aspx?ID=6664&productid=2723>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>