

## Purified anti-mouse IL-6

Rat IgG1, kappa  
MP5-20F3



**BioLegend**  
From Revvity

500 µg  
Artikel-Nr.: B504502 | **BioLegend** | Hersteller-Nr.: 504502

**111,00 €\***

\*zzgl. MwSt. zzgl. Versandkosten

### Beschreibung

**App\_Abbrev:** FA, FC, ICC, ICFC, IHC, IP, WB  
**Clone:** MP5-20F3  
**Format:** PURE  
**Isotype:** Rat IgG1, kappa  
**Species of Origin:** Rat  
**Verpackung:** 500 µg

Purified anti-mouse IL-6

Clone: MP5-20F3

Isotype: Rat IgG1, kappa

Other\_Names: Interleukin-6, B cell stimulating factor-2 (BSF-2), Cytotoxic T cell differentiation factor (CDF), Hepatocyte stimulating factor (HSF), Hybridoma/plasmacytoma growth factor (HPGF)

Size: 500 µg

Reactivity: Mouse

Format: PURE

App\_Abbrev: FC, ICFC, ICC, IHC, IP, WB, FA

Species of Origin: Rat

Weitere Informationen zum Produkt. Link zum Hersteller: <https://www.biolegend.com/Default.aspx?ID=6664&productid=973>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>