

PE anti-mouse IL-6

Rat IgG1, kappa
MP5-20F3



BioLegend
From Revvity

25 µg
Artikel-Nr.: B504503 | **BioLegend** | Hersteller-Nr.: 504503

77,00 €*

*zzgl. MwSt. zzgl. Versandkosten

Beschreibung

App_Abbrev: FC, ICFC
Clone: MP5-20F3
Format: PE
Isotype: Rat IgG1, kappa
Species of Origin: Rat
Verpackung: 25 µg

PE anti-mouse IL-6

Clone: MP5-20F3

Isotype: Rat IgG1, kappa

Other_Names: Interleukin-6, B cell stimulating factor-2 (BSF-2), Cytotoxic T cell differentiation factor (CDF), Hepatocyte stimulating factor (HSF), Hybridoma/plasmacytoma growth factor (HPGF)

Size: 25 µg

Reactivity: Mouse

Format: PE

App_Abbrev: FC, ICFC

Species of Origin: Rat

Weitere Informationen zum Produkt. Link zum Hersteller: <https://www.biolegend.com/Default.aspx?ID=6664&productid=972>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>