

Purified anti-mouse IL-6

Rat IgG2a, kappa
MP5-32C11



BioLegend
From Revvity

500 µg
Artikel-Nr.: B504604 | **BioLegend** | Hersteller-Nr.: 504604

188,00 €*

*zzgl. MwSt. zzgl. Versandkosten

Beschreibung

App_Abbrev: FA, FC, ICC, ICFC, IHC, IP, WB
Clone: MP5-32C11
Format: PURE
Isotype: Rat IgG2a, kappa
Species of Origin: Rat
Verpackung: 500 µg

Purified anti-mouse IL-6

Clone: MP5-32C11

Isotype: Rat IgG2a, kappa

Other_Names: Interleukin-6, B cell stimulating factor-2 (BSF-2), Cytotoxic T cell differentiation factor (CDF), Hepatocyte stimulating factor (HSF), Hybridoma/plasmacytoma growth factor (HPGF)

Size: 500 µg

Reactivity: Mouse

Format: PURE

App_Abbrev: FC, ICFC, ICC, IHC, IP, WB, FA

Species of Origin: Rat

Weitere Informationen zum Produkt. Link zum Hersteller:

<https://www.biolegend.com/Default.aspx?ID=6664&productid=3836>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>