

Purified anti-mouse IL-9

Rat IgG2a, kappa
D8402E8



BioLegend
From Revvity

500 µg
Artikel-Nr.: B504702 | **BioLegend** | Hersteller-Nr.: 504702

198,00 €*

*zzgl. MwSt. zzgl. Versandkosten

Beschreibung

App_Abbrev: FA, FC, ICC, ICFC, IHC, IP, WB
Clone: D8402E8
Format: PURE
Isotype: Rat IgG2a, kappa
Species of Origin: Rat
Verpackung: 500 µg

Purified anti-mouse IL-9

Clone: D8402E8
Isotype: Rat IgG2a, kappa
Other_Names: Interleukin-9, P40, Mast cell growth factor (MCGF), Megakaryoblast growth factor, T cell growth factor-3 (TCGF-3)
Size: 500 µg

Reactivity: Mouse
Format: PURE
App_Abbrev: FC, ICFC, ICC, IHC, IP, WB, FA
Species of Origin: Rat

Weitere Informationen zum Produkt. Link zum Hersteller: <https://www.biolegend.com/Default.aspx?ID=6664&productid=932>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>