

Purified anti-mouse IL-10

Rat IgG1, kappa
JES5-2A5



BioLegend
From Revvity

50 µg
Artikel-Nr.: B504901 | [BioLegend](#) | Hersteller-Nr.: 504901

67,00 €*

*zzgl. MwSt. [zzgl. Versandkosten](#)

Beschreibung

App_Abbrev: FA, FC, ICC, ICFC, IHC, IP, WB
Clone: JES5-2A5
Format: PURE
Isotype: Rat IgG1, kappa
Species of Origin: Rat
Verpackung: 50 µg

Purified anti-mouse IL-10

Clone: JES5-2A5

Isotype: Rat IgG1, kappa

Other_Names: Interleukin-10, Cytokine synthesis inhibitory factor (CSIF), B cell derived T cell growth factor (B-TCGF), T cell growth inhibitory factor (TGIF)

Size: 50 µg

Reactivity: Mouse

Format: PURE

App_Abbrev: FC, ICFC, ICC, IHC, IP, WB, FA

Species of Origin: Rat

Weitere Informationen zum Produkt. Link zum Hersteller: <https://www.biolegend.com/Default.aspx?ID=6664&productid=947>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>