

## PE anti-mouse IL-10

Rat IgG2b, kappa  
JES5-16E3



**BioLegend**  
From Revvity

25 µg  
Artikel-Nr.: B505007 | **BioLegend** | Hersteller-Nr.: 505007

**68,00 €\***

\*zzgl. MwSt. zzgl. Versandkosten

### Beschreibung

**App\_Abbrev:** FC, ICFC  
**Clone:** JES5-16E3  
**Format:** PE  
**Isotype:** Rat IgG2b, kappa  
**Species of Origin:** Rat  
**Verpackung:** 25 µg

PE anti-mouse IL-10

Clone: JES5-16E3  
Isotype: Rat IgG2b, kappa  
Other\_Names: Interleukin-10, Cytokine synthesis inhibitory factor (CSIF), B cell derived T cell growth factor (B-TCGF)  
Size: 25 µg

Reactivity: Mouse  
Format: PE  
App\_Abbrev: FC, ICFC  
Species of Origin: Rat

Weitere Informationen zum Produkt. Link zum Hersteller: <https://www.biolegend.com/Default.aspx?ID=6664&productid=944>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>