

APC anti-mouse IL-10

Rat IgG2b, kappa
JES5-16E3



BioLegend
From Revvity

25 µg
Artikel-Nr.: B505009 | [BioLegend](#) | Hersteller-Nr.: 505009

82,00 €*

*zzgl. MwSt. [zzgl. Versandkosten](#)

Beschreibung

App_Abbrev: FC, ICFC
Clone: JES5-16E3
Format: APC
Isotype: Rat IgG2b, kappa
Species of Origin: Rat
Verpackung: 25 µg

APC anti-mouse IL-10

Clone: JES5-16E3

Isotype: Rat IgG2b, kappa

Other_Names: Interleukin-10, Cytokine synthesis inhibitory factor (CSIF), B cell derived T cell growth factor (B-TCGF)

Size: 25 µg

Reactivity: Mouse

Format: APC

App_Abbrev: FC, ICFC

Species of Origin: Rat

Weitere Informationen zum Produkt. Link zum Hersteller: <https://www.biolegend.com/Default.aspx?ID=6664&productid=941>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>