

Purified anti-mouse IL-12/IL-23 p40

Rat IgG1, kappa
C15.6



BioLegend
From Revvity

50 µg
Artikel-Nr.: B505201 | **BioLegend** | Hersteller-Nr.: 505201

67,00 €*

*zzgl. MwSt. zzgl. Versandkosten

Beschreibung

App_Abbrev: FA, FC, ICC, ICFC, IHC, IP, WB
Clone: C15.6
Format: PURE
Isotype: Rat IgG1, kappa
Species of Origin: Rat
Verpackung: 50 µg

Purified anti-mouse IL-12/IL-23 p40

Clone: C15.6

Isotype: Rat IgG1, kappa

Other_Names: Interleukin-12 p40, Interleukin-23 p40, Cytotoxic lymphocyte maturation factor (CLMF), Natural killer cell stimulatory factor (NKSF), CTL maturation factor (TcMF), T-cell stimulating factor (TSF)

Size: 50 µg

Reactivity: Mouse

Format: PURE

App_Abbrev: FC, ICFC, ICC, IHC, IP, WB, FA

Species of Origin: Rat

Weitere Informationen zum Produkt. Link zum Hersteller: <https://www.biolegend.com/Default.aspx?ID=6664&productid=929>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>