

APC anti-mouse IL-12/IL-23 p40

Rat IgG1, kappa
C15.6



BioLegend
From Revvity

100 µg
Artikel-Nr.: B505206 | [BioLegend](#) | Hersteller-Nr.: 505206

263,00 €*

*zzgl. MwSt. [zzgl. Versandkosten](#)

Beschreibung

App_Abbrev: FC, ICFC
Clone: C15.6
Format: APC
Isotype: Rat IgG1, kappa
Species of Origin: Rat
Verpackung: 100 µg

APC anti-mouse IL-12/IL-23 p40

Clone: C15.6

Isotype: Rat IgG1, kappa

Other_Names: Interleukin-12 p40, Interleukin-23 p40, Cytotoxic lymphocyte maturation factor (CLMF), Natural killer cell stimulatory factor (NKSF), CTL maturation factor (TcMF), T-cell stimulating factor (TSF)

Size: 100 µg

Reactivity: Mouse

Format: APC

App_Abbrev: FC, ICFC

Species of Origin: Rat

Weitere Informationen zum Produkt. Link zum Hersteller: <https://www.biolegend.com/Default.aspx?ID=6664&productid=926>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>