

## Purified anti-mouse GM-CSF

Rat IgG2a, kappa  
MP1-22E9



**BioLegend**  
From Revvity

500 µg  
Artikel-Nr.: B505402 | **BioLegend** | Hersteller-Nr.: 505402

**111,00 €\***

\*zzgl. MwSt. zzgl. Versandkosten

### Beschreibung

**App\_Abbrev:** ELISA Capture, ELISPOT Capture, ICC, IHC-F, Neut  
**Clone:** MP1-22E9  
**Format:** PURE  
**Isotype:** Rat IgG2a, kappa  
**Species of Origin:** Rat  
**Verpackung:** 500 µg

Purified anti-mouse GM-CSF

Clone: MP1-22E9

Isotype: Rat IgG2a, kappa

Other\_Names: Granulocyte/macrophage-colony stimulating factor, CSF-alpha, Pluripoietin-alpha, Eosinophil colony stimulating factor (Eo-CSF), Burst promoting activity (BPA)

Size: 500 µg

Reactivity: Mouse

Format: PURE

App\_Abbrev: ELISA Capture, ELISPOT Capture, Neut, ICC, IHC-F

Species of Origin: Rat

Weitere Informationen zum Produkt. Link zum Hersteller: <https://www.biolegend.com/Default.aspx?ID=6664&productid=959>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>