

Purified anti-mouse/rat IL-1beta

Armenian Hamster IgG
B122



BioLegend
From Revvity

500 µg
Artikel-Nr.: B503502 | **BioLegend** | Hersteller-Nr.: 503502

115,00 €*

*zzgl. MwSt. zzgl. Versandkosten

Beschreibung

App_Abbrev: FA, FC, ICC, ICFC, IHC, IP, WB
Clone: B122
Format: PURE
Isotype: Armenian Hamster IgG
Species of Origin: Armenian Hamster
Verpackung: 500 µg

Purified anti-mouse/rat IL-1beta

Clone: B122

Isotype: Armenian Hamster IgG

Other_Names: Interleukin-1beta, Catabolin, Hematopoietin-1 (H1), IFNbeta-inducing factor, Interleukin-beta, Osteoclast activating factor (OAF)

Size: 500 µg

Reactivity: Mouse, Rat

Format: PURE

App_Abbrev: FC, ICFC, ICC, IHC, IP, WB, FA

Species of Origin: Armenian Hamster

Weitere Informationen zum Produkt. Link zum Hersteller: <https://www.biolegend.com/Default.aspx?ID=6664&productid=908>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>