

Purified anti-mouse IL-2

Rat IgG2a, kappa
JES6-1A12



BioLegend
From Revvity

50 µg
Artikel-Nr.: B503701 | [BioLegend](#) | Hersteller-Nr.: 503701

67,00 €*

*zzgl. MwSt. [zzgl. Versandkosten](#)

Beschreibung

App_Abbrev: FA, FC, ICC, ICFC, IHC, IP, WB
Clone: JES6-1A12
Format: PURE
Isotype: Rat IgG2a, kappa
Species of Origin: Rat
Verpackung: 50 µg

Purified anti-mouse IL-2

Clone: JES6-1A12

Isotype: Rat IgG2a, kappa

Other_Names: Interleukin-2, T cell growth factor (TCGF), Eosinophil differentiation factor (EDF), Killer cell helper factor (KHF), Macrophage-activating factor for cytotoxicity I (MAF-C I), Thymocyte differentiation factor (TDF)

Size: 50 µg

Reactivity: Mouse

Format: PURE

App_Abbrev: FC, ICFC, ICC, IHC, IP, WB, FA

Species of Origin: Rat

Weitere Informationen zum Produkt. Link zum Hersteller: <https://www.biolegend.com/Default.aspx?ID=6664&productid=949>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>