

Alexa Fluor® 647 anti-mouse IL-2

Rat IgG2b, kappa
JES6-5H4



BioLegend
From Revvity

100 µg

Artikel-Nr.: B503814 | [BioLegend](#) | Hersteller-Nr.: 503814

212,00 €*

*zzgl. MwSt. [zzgl. Versandkosten](#)

Beschreibung

App_Abbrev: FC, ICFC
Clone: JES6-5H4
Format: A647
Isotype: Rat IgG2b, kappa
Species of Origin: Rat
Verpackung: 100 µg

Alexa Fluor® 647 anti-mouse IL-2

Clone: JES6-5H4

Isotype: Rat IgG2b, kappa

Other_Names: Interleukin-2, T cell growth factor (TCGF), Eosinophil differentiation factor (EDF), Killer cell helper factor (KHF), Macrophage-activating factor for cytotoxicity I (MAF-C I), Thymocyte differentiation factor (TDF)

Size: 100 µg

Reactivity: Mouse

Format: A647

App_Abbrev: FC, ICFC

Species of Origin: Rat

Weitere Informationen zum Produkt. Link zum Hersteller:

<https://www.biolegend.com/Default.aspx?ID=6664&productid=2716>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>