

Purified anti-mouse IL-3

Rat IgG1, kappa
MP2-8F8



BioLegend
From Revvity

500 µg
Artikel-Nr.: B503902 | [BioLegend](#) | Hersteller-Nr.: 503902

111,00 €*

*zzgl. MwSt. [zzgl. Versandkosten](#)

Beschreibung

App_Abbrev: FA, FC, ICC, ICFC, IHC, IP, WB
Clone: MP2-8F8
Format: PURE
Isotype: Rat IgG1, kappa
Species of Origin: Rat
Verpackung: 500 µg

Purified anti-mouse IL-3

Clone: MP2-8F8

Isotype: Rat IgG1, kappa

Other_Names: Interleukin-3, Burst promoting activity, Eosinophil colony stimulating factor (Eo-CSF), Hematopoietic cell growth factor (HCGF), Mast (MGF/MCGF), Multi-colony stimulating (Multi-CSF), P cell stimulating activity (PCSA), Thy1 inducing factor

Size: 500 µg

Reactivity: Mouse

Format: PURE

App_Abbrev: FC, ICFC, ICC, IHC, IP, WB, FA

Species of Origin: Rat

Weitere Informationen zum Produkt. Link zum Hersteller: <https://www.biolegend.com/Default.aspx?ID=6664&productid=964>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>