

PE anti-mouse IL-3

Rat IgG1, kappa
MP2-8F8



BioLegend
From Revvity

25 µg
Artikel-Nr.: B503903 | **BioLegend** | Hersteller-Nr.: 503903

88,00 €*

*zzgl. MwSt. zzgl. Versandkosten

Beschreibung

App_Abbrev: FC, ICFC
Clone: MP2-8F8
Format: PE
Isotype: Rat IgG1, kappa
Species of Origin: Rat
Verpackung: 25 µg

PE anti-mouse IL-3

Clone: MP2-8F8

Isotype: Rat IgG1, kappa

Other_Names: Interleukin-3, Burst promoting activity, Eosinophil colony stimulating factor (Eo-CSF), Hematopoietic cell growth factor (HCGF), Mast (MGF/MCGF), Multi-colony stimulating (Multi-CSF), P cell stimulating activity (PCSA), Thy1 inducing factor

Size: 25 µg

Reactivity: Mouse

Format: PE

App_Abbrev: FC, ICFC

Species of Origin: Rat

Weitere Informationen zum Produkt. Link zum Hersteller: <https://www.biolegend.com/Default.aspx?ID=6664&productid=963>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>