

Purified anti-mouse IL-4

Rat IgG1, kappa
11B11



BioLegend
From Revvity

50 µg
Artikel-Nr.: B504101 | [BioLegend](#) | Hersteller-Nr.: 504101

67,00 €*

*zzgl. MwSt. [zzgl. Versandkosten](#)

Beschreibung

App_Abbrev: FA, FC, ICC, ICFC, IHC, IP, WB
Clone: 11B11
Format: PURE
Isotype: Rat IgG1, kappa
Species of Origin: Rat
Verpackung: 50 µg

Purified anti-mouse IL-4

Clone: 11B11

Isotype: Rat IgG1, kappa

Other_Names: Interleukin-4, Ia inducing factor (IaIF), B cell stimulating factor-1 (BSF-1), Hodgkin's cell growth factor (HCGF), Mast cell growth factor-2 (MCGF-2), Macrophage fusion factor (MFF), T cell growth factor-2 (TCGF-2)

Size: 50 µg

Reactivity: Mouse

Format: PURE

App_Abbrev: FC, ICFC, ICC, IHC, IP, WB, FA

Species of Origin: Rat

Weitere Informationen zum Produkt. Link zum Hersteller: <https://www.biolegend.com/Default.aspx?ID=6664&productid=894>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>