

Alexa Fluor® 488 anti-mouse IL-4

Rat IgG1, kappa
11B11



BioLegend
From Revvity

100 µg

Artikel-Nr.: B504109 | BioLegend | Hersteller-Nr.: 504109

215,00 €*

*zzgl. MwSt. zzgl. Versandkosten

Beschreibung

App_Abbrev: FC, ICFC
Clone: 11B11
Format: A488
Isotype: Rat IgG1, kappa
Species of Origin: Rat
Verpackung: 100 µg

Alexa Fluor® 488 anti-mouse IL-4

Clone: 11B11

Isotype: Rat IgG1, kappa

Other_Names: Interleukin-4, Ia inducing factor (IaIF), B cell stimulating factor-1 (BSF-1), Hodgkin's cell growth factor (HCGF), Mast cell growth factor-2 (MCGF-2), Macrophage fusion factor (MFF), T cell growth factor-2 (TCGF-2)

Size: 100 µg

Reactivity: Mouse

Format: A488

App_Abbrev: FC, ICFC

Species of Origin: Rat

Weitere Informationen zum Produkt. Link zum Hersteller:

<https://www.biolegend.com/Default.aspx?ID=6664&productid=2717>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>