

## Purified anti-mouse/human IL-5

Rat IgG1, kappa  
TRFK5



**BioLegend**  
From Revvity

500 µg  
Artikel-Nr.: B504302 | [BioLegend](#) | Hersteller-Nr.: 504302

**111,00 €\***

\*zzgl. MwSt. [zzgl. Versandkosten](#)

### Beschreibung

**App\_Abbrev:** FA, FC, ICC, ICFC, IHC, IP, WB  
**Clone:** TRFK5  
**Format:** PURE  
**Isotype:** Rat IgG1, kappa  
**Species of Origin:** Rat  
**Verpackung:** 500 µg

Purified anti-mouse/human IL-5

Clone: TRFK5

Isotype: Rat IgG1, kappa

Other\_Names: Interleukin-5, Eosinophil differentiation factor (EDF), Eosinophil colony stimulating factor (Eo-CSF), B cell growth factor-2 (BCGF-2), B cell differentiation factor for IgM (BCDF-m), IgA enhancing factor (IgA-EF), T-cell replacing factor-1 (TRF-1)

Size: 500 µg

Reactivity: Mouse, Human

Format: PURE

App\_Abbrev: FC, ICFC, ICC, IHC, IP, WB, FA

Species of Origin: Rat

Weitere Informationen zum Produkt. Link zum Hersteller: <https://www.biolegend.com/Default.aspx?ID=6664&productid=992>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>