

PE anti-mouse/human IL-5

Rat IgG1, kappa
TRFK5



BioLegend
From Revvity

25 µg
Artikel-Nr.: B504303 | [BioLegend](#) | Hersteller-Nr.: 504303

77,00 €*

*zzgl. MwSt. [zzgl. Versandkosten](#)

Beschreibung

App_Abbrev: FC, ICFC
Clone: TRFK5
Format: PE
Isotype: Rat IgG1, kappa
Species of Origin: Rat
Verpackung: 25 µg

PE anti-mouse/human IL-5

Clone: TRFK5

Isotype: Rat IgG1, kappa

Other_Names: Interleukin-5, Eosinophil differentiation factor (EDF), Eosinophil colony stimulating factor (Eo-CSF), B cell growth factor-2 (BCGF-2), B cell differentiation factor for IgM (BCDF-m), IgA enhancing factor (IgA-EF), T-cell replacing factor-1 (TRF-1)

Size: 25 µg

Reactivity: Mouse, Human

Format: PE

App_Abbrev: FC, ICFC

Species of Origin: Rat

Weitere Informationen zum Produkt. Link zum Hersteller: <https://www.biolegend.com/Default.aspx?ID=6664&productid=991>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>