

## PE anti-mouse IL-1alpha

Armenian Hamster IgG  
ALF-161



**BioLegend**  
From Revvity

25 µg  
Artikel-Nr.: B503203 | **BioLegend** | Hersteller-Nr.: 503203

**88,00 €\***

\*zzgl. MwSt. zzgl. Versandkosten

### Beschreibung

**App\_Abbrev:** FC, ICFC  
**Clone:** ALF-161  
**Format:** PE  
**Isotype:** Armenian Hamster IgG  
**Species of Origin:** Armenian Hamster  
**Verpackung:** 25 µg

PE anti-mouse IL-1alpha

Clone: ALF-161

Isotype: Armenian Hamster IgG

Other\_Names: Interleukin-1alpha, Lymphocyte-activating factor (LAF), Endogenous pyrogen (EP), Leukocyte endogenous mediator (LEM), Mononuclear cell factor (MCF/MNCF)

Size: 25 µg

Reactivity: Mouse

Format: PE

App\_Abbrev: FC, ICFC

Species of Origin: Armenian Hamster

Weitere Informationen zum Produkt. Link zum Hersteller: <https://www.biolegend.com/Default.aspx?ID=6664&productid=903>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>