

Purified anti-human LT-alpha (TNF-beta)

Mouse IgG1, kappa
359-238-8



BioLegend
From Revvity

500 µg
Artikel-Nr.: B503002 | **BioLegend** | Hersteller-Nr.: 503002

188,00 €*

*zzgl. MwSt. zzgl. Versandkosten

Beschreibung

App_Abbrev: FA, FC, ICC, ICFC, IHC, IP, WB
Clone: 359-238-8
Format: PURE
Isotype: Mouse IgG1, kappa
Species of Origin: Mouse
Verpackung: 500 µg

Purified anti-human LT-alpha (TNF-beta)

Clone: 359-238-8

Isotype: Mouse IgG1, kappa

Other_Names: Lymphotoxin-alpha (LT- alpha), Tumor necrosis factor-beta (TNF-beta), Coley's toxin, Hemorrhagic factor, Necrosin, Natural killer cytotoxic factor (NKCF), Differentiation inducing factor (DIF), TNFSF-1

Size: 500 µg

Reactivity: Human

Format: PURE

App_Abbrev: FC, ICFC, ICC, IHC, IP, WB, FA

Species of Origin: Mouse

Weitere Informationen zum Produkt. Link zum Hersteller:

<https://www.biolegend.com/Default.aspx?ID=6664&productid=1319>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>