

Purified anti-mouse IL-1alpha

Armenian Hamster IgG
ALF-161



BioLegend
From Revvity

500 µg
Artikel-Nr.: B503202 | BioLegend | Hersteller-Nr.: 503202

109,00 €*

*zzgl. MwSt. [zzgl. Versandkosten](#)

Beschreibung

App_Abbrev: FA, FC, ICC, ICFC, IHC, IP, WB
Clone: ALF-161
Format: PURE
Isotype: Armenian Hamster IgG
Species of Origin: Armenian Hamster
Verpackung: 500 µg

Purified anti-mouse IL-1alpha

Clone: ALF-161

Isotype: Armenian Hamster IgG

Other_Names: Interleukin-1alpha, Lymphocyte-activating factor (LAF), Endogenous pyrogen (EP), Leukocyte endogenous mediator (LEM), Mononuclear cell factor (MCF/MNCF)

Size: 500 µg

Reactivity: Mouse

Format: PURE

App_Abbrev: FC, ICFC, ICC, IHC, IP, WB, FA

Species of Origin: Armenian Hamster

Weitere Informationen zum Produkt. Link zum Hersteller: <https://www.biolegend.com/Default.aspx?ID=6664&productid=904>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>