

Purified anti-human IL-10

Rat IgG2a, kappa
JES3-12G8



BioLegend
From Revvity

500 µg

Artikel-Nr.: B501504 | [BioLegend](#) | Hersteller-Nr.: 501504

188,00 €*

*zzgl. MwSt. [zzgl. Versandkosten](#)

Beschreibung

App_Abbrev: FA, FC, ICC, ICFC, IHC, IP, WB
Clone: JES3-12G8
Format: PURE
Isotype: Rat IgG2a, kappa
Species of Origin: Rat
Verpackung: 500 µg

Purified anti-human IL-10

Clone: JES3-12G8

Isotype: Rat IgG2a, kappa

Other_Names: Interleukin-10, B cell derived T cell growth factor (B-TCGF), Cytokine synthesis inhibitory factor (CSIF), T-cell growth inhibitory factor (TGIF)

Size: 500 µg

Reactivity: Human

Format: PURE

App_Abbrev: FC, ICFC, ICC, IHC, IP, WB, FA

Species of Origin: Rat

Weitere Informationen zum Produkt. Link zum Hersteller:

<https://www.biolegend.com/Default.aspx?ID=6664&productid=4560>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>